








## DNA Extraction



**Goal: Visualize the physical nature of DNA**

### Materials:

-  • .9% NaCl (saline)
-  • Dixie cups
-  • 10 ml graduated cylinder
-  • Test tube and rack
-  • 25% soapy water solution
-  • Ethanol (in a squeeze bottle)
- Toothpicks
- Strawberries and Kiwis
- Ziploc baggies
- Table salt (dry NaCl)
- Gauze to filter smushed fruit mix
-  • Funnels
- Scissors



**Step One – Cheek Cell DNA Extraction (*Complete this stage before moving on!*) Procedure:**

- Vigorously swish 10ml of .9% saline solution (9 grams of dry NaCl in 1 liter of water) in your mouth for 30-60 seconds. The Dixie cups have 10ml of saline solution in them already.


 Think & Write about this: Why did you use saline instead of plain water? 

- Gently spit the saline/spit mix into your cup.
- Select a test tube. Note the number \_\_\_\_\_. You'll use this test tube again.
- Carefully pour the saline/spit mix into your clean test tube and place it in a test tube rack.
- Measure 5 ml of soapy water and carefully pour into the test tube.

→ Hint: use the graduated cylinder that's marked for soap so you don't have to wash it.

 Think & Write about this: Why did you use soap? 

- Place your thumb over the test tube opening and gently tilt the test tube back and forth several times (10 times works well) to mix well.
- Using a squeeze bottle, run about 1/2" of chilled ethanol down the inside edge of the test tube – be careful so the two liquids don't mix. You're trying to create two separate layers.
- You should observe (see) a white stringy substance forming at the place where the two layers meet. Don't be alarmed if some of the material rises to the top of the ethanol – that's okay.
- Twirl a toothpick or wooden stirrer in the DNA to extract the DNA. The "strings" will wind around the toothpick/stirrer. Cool, hunh?

 Think & Write about this: What are you seeing? What does it look like?


- Wash your test tube and take it with you to the next station (there's more work to be done!).

### Step Two – Fruit DNA Extraction (using Strawberries and Kiwi) Procedure:

- Select a piece of fruit (always a good place to start).


Hint: Strawberries work really well because most of them are octoploid. This means that they have eight (*octo means eight*) copies of each chromosome. In comparison, humans are diploid (*di means two*), having two copies of each chromosome.

- Place a small (1-inch cube) piece of the chosen fruit into a Ziploc baggie Press the air out and seal the bag. Mash the bagged fruit **really well** with your fingers.

 Think & Write about this: Why are you mashing the fruit?



- Add 10ml of the 25 % soap solution and a pinch of salt to the baggie. Press the air out and seal the bag. Mash the fruit, soap and salt for another minute.
- Place a gauze or filter-lined funnel into your test tube while it's still in the test tube rack.
- Use scissors to cut a corner of the baggie and squeeze the fruit/soap/salt mixture into a gauze or filter-lined funnel, collecting the liquid in your test tube.
- Throw away the fruit pulp and filter in the trash or provided container.
- Slowly drip about a 1/2" of chilled ethanol along the side of the test tube using the squeeze bottle. Remember: the goal is to form a new layer of ethanol on top of the liquid.
- You should observe a white filamentous (stringy) substance forming at the place where the two layers meet. Did you get more DNA material this time? Why/why not, do you think?
- Twirl a toothpick/wooden stirrer in the DNA to extract the DNA. The "strings" will wind around the toothpick. Still cool, hunh?
- Wash your test tube!

 Think & Write about this: Did your cheek cell DNA extraction look like the fruit DNA? Why do you think it looked the same and/or different?